



Triangle:  
Sign of PIN(1)



- |    |                          |
|----|--------------------------|
| 1  | PIN(1): Out GND          |
| 2  | PIN(2): Output 1, or 9   |
| 3  | PIN(3): Out GND          |
| 4  | PIN(4): Output 2, or 10  |
| 5  | PIN(5): Out GND          |
| 6  | PIN(6): Output 3, or 11  |
| 7  | PIN(7): Out GND          |
| 8  | PIN(8): Output 4, or 12  |
| 9  | PIN(9): Out GND          |
| 10 | PIN(10): Output 5, or 13 |
| 11 | PIN(11): Out GND         |
| 12 | PIN(12): Output 6, or 14 |
| 13 | PIN(13): Out GND         |
| 14 | PIN(14): Output 7, or 15 |
| 15 | PIN(15): Out GND         |
| 16 | PIN(16): Output 8, or 16 |
| 17 | PIN(17): Out GND         |
| 18 | PIN(18): Out GND         |
| 19 | PIN(19): Out GND         |
| 20 | PIN(20): Out GND         |

The appropriate 20-pin, 2.54 mm female plug compatible with this output connector is IDS-020-S100-01/P manufactured by E-tec Corp.

## Output Connector of each 8-channel sections of BioAmp, V.1